



lsxisxray@aol.com

1-800-777-4518

### Veterinary 300ma and 100ma Technique Chart

In-Table Procedure																
Dog Thorax 300 ma 1/60 sec	CM	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
	KVP	64	66	68	70	72	74	76	78	80	82	84	86	88	92	94
Dog Abdomen 300 ma 1/30 sec	CM	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
	KVP	62	66	68	70	72	74	76	80	82	84	86	88	94	96	

Table-Top Procedure																
Cat Abdomen 300 ma 1/30 sec	CM	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
	KVP					50	54	58	60	62	62	64	66	68	70	74
Head/Spine/Pelvis/Should 300 ma 1/20 sec	CM	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
	KVP	64	66	68	70	70	72	74	76	78	80	84	86			

In-Table & Table-Top Procedure												
Cat Thorax w/Grid 300 ma 1/60 sec	CM	2	3	4	5	6	7	8	9	10	11	12
	KVP				48	50	52	54	56	58	62	64
Extremity & Skull (Table-T 100 ma 1/40 sec	CM	2	3	4	5	6	7	8	9	10	11	12
	KVP	52	54	56	58	60	62	64	66	68	70	72



porated

g

*Quality*

1-800-777-4518

25
98
25

16
76
25

