

## Veterinary Technique Chart

Cassette/Film - 400 Speed System

### In-Table Procedure

Dog Thorax	CM	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
300 ma 1/60 sec	kvp	64	66	68	70	72	74	76	78	80	82	84	86	88	92	94	98
Dog Abdomen	CM	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
300 ma 1/30 sec	kvp	62	66	68	70	72	74	76	80	82	84	86	88	94	96		

### Table-Top Procedure

Cat Abdomen	CM	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
300 ma 1/30 sec	kvp					50	54	58	60	62	62	64	66	68	70	74	76
Head/Spine/Pelvis/Shoulder	CM	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
300 ma 1/20 sec	kvp	64	66	68	70	70	72	74	76	78	80	84	86				

### In-Table & Table-Top Procedure

Cat Thorax w/Grid	CM	2	3	4	5	6	7	8	9	10	11	12					
300 ma 1/60 sec	kvp				48	50	52	54	56	58	62	64					
Extremity & Skull (Table-Top)	CM	2	3	4	5	6	7	8	9	10	11	12					
100 ma 1/40 sec	kvp	52	54	56	58	60	62	64	66	68	70	72					